



Manage Your “Sugar” with An Easy Touch

Limpid

Blood Glucose Meter

Blood glucose testing for professional
and self-use

Small size with large display

0.5 μ L blood taken, least pain

Accurate result within 5 seconds



Product Specifications

Parameter	Specification
Test Items	Blood glucose
Sample Type	Capillary whole blood
Sample Volume	0.5 μ L
Testing Time	5 seconds
Display unit	mmol/L or mg/dL
Test Range	0.6~33.3mmol/L(10~600mg/dL)
Hematocrit range	15~70%
Strip Storage Condition	Temperature:2~30°C; humidity: $\leq$ 90%;
Meter Size	86.9mm×53.4mm×18.9mm
Meter Weight	49g (without batteries)
Analyzer Power Supply	CR2032 batteries (3V)
Working Condition	Temperature: 5~45 °C ; Humidity: $\leq$ 90%
Data Storage	900
Date Transfer and Connectivity	Micro USB



Limpid
Blood Glucose Meter

Order Information

Product Name	Catalog No.	Components
Limpid Blood Glucose Meter	GLU-M11	1 Blood glucose meter 2 CR2032 batteries (3V) 1 Lancing Device 1 Carry Case 1 User Manual 1 Quick Start Guide 1 Warranty Card

Product Name	Catalog No.	Package
Limpid Blood Glucose Test Strip	GLU-S42	50 Strips/Kit