


Home Test

Visual Cholesterol Test Kit

• QUCARE Self Total Cholesterol

Simple and easy-to-use test kit without any instrument. It is a screening test to provide a total cholesterol level for everyone in 2 minutes at anytime and anywhere. Keep a healthy lifestyle with a quick test.




Easy use
Anytime Anywhere
(Not required to send sample to Lab)



Quick test
in 2 minutes



Simple way
without instrument




How to use a test kit

Step 1)



- Wash and Dry your hands

Step 2)



- Prick your finger and get a blood sample (only 10uL)

Step 3)



- Apply the blood sample and check the test result after 2 minutes.

Technical Specification

- 99-250 mg/dl(2.6 ~6.5 mmol/l)
- Capillary whole blood
- 2 minutes
- 10 uL(Capillary included)

